


C2F Components

1. C2F Assembly (required for first case; 1 assembly per implant site)
Includes the following:
 - a. C2F Blockout Rod
 - b. C2F Analog
 - c. C2F MUA
 - d. C2F Temporary Coping
 - e. C2F Gold Screw
 - f. C2F Drill
2. C2F Model(s)
3. C2F Provisional(s)



Conversion Tools

1. Straight handpiece with foot pedal
2. Acrylic resin (ex. STELLAR) with multiple applicator tips
3. RAPID Appliance
4. Superglue and accelerant
5. Vertex Articulator
6. Straight small hex driver
7. Straight prosthetic screwdriver
8. CHROME Conversion Bur Kit or similar kit with acrylic adjustment burs and polishing wheels
9. Curing light if using STELLAR DC Acrylic Dual Cure



To learn more, please visit
www.roedentallab.com/c2f